

## DATASHEET

ZEOtube TA 55mm 7" 100 PC



Glass Tube for NMR analysis

TA - Throw Away

Release date: 18.04.2023

Replaces all other versions

Description	
Material	Borosilicate Glass 5.1 ASTM Type 1 class B glass
Article Number	400490
Tubes per Pack	100 tubes with orange caps
Outer Diameter	$4.950 \pm 0.030$ mm
Wall thickness	0.43 mm
Concentricity	25 $\mu$ m
Camber	75 $\mu$ m
Length	7 inch

ZEOtube

