

DATASHEET

ZEOtube HP 5mm 7" 5 PC



Glass Tube for NMR analysis

HP - High Precision

Release date: 18.04.2023

Replaces all other versions

Description	
Material	Borosilicate Glass 5.1 ASTM Type 1 class B glass
Article Number	400491
Tubes per Pack	5 tubes with sky blue caps
Outer Diameter	4.950 ± 0.015 mm
Wall thickness	0.43 mm
Concentricity	7 μ m
Camber	13 μ m
Length	7 inch

ZEOtube

