

## DATASHEET

ZEOtube UP 5mm 7" 5 PC



Glass Tube for NMR analysis  
UP - Ultra Precision

Release date: 18.04.2023

Replaces all other versions

Description	
Material	Borosilicate Glass 5.1 ASTM Type 1 class B glass
Article Number	400492
Tubes per Pack	5 tubes with black caps
Outer Diameter	$4.950 \pm 0.010$ mm
Wall thickness	0.43 mm
Concentricity	4 $\mu$ m
Camber	4 $\mu$ m
Length	7 inch

